



## Ziwig Endotest<sup>®</sup>

A saliva test for the early and reliable diagnosis  
of endometriosis

Women suffering from endometriosis, for the most part, experience a long delay between the onset of symptoms and diagnosis (between 7 and 11 years on average)<sup>1,2</sup>. During this time they may see multiple physicians, undergo numerous examinations and, in some cases, have surgery before a diagnosis is reached.

A non-invasive test now allows women suffering from symptoms of endometriosis to avoid this long and tortuous medical journey. A team, made up of expert doctors specialized in endometriosis, engineers in artificial intelligence and medical researchers has recently developed a saliva based diagnostic test that uses Next Generation Sequencing of microRNAs.



This simple and non-invasive test, called Ziwig Endotest®, allows for the early diagnosis of all types of endometriosis with a performance that exceeds that of conventional diagnostic tests<sup>3-8</sup> (sensitivity > 97%, specificity > 93%).

Ziwig Endotest® is based on robust clinical evidence including validation by one of the largest clinical trials ever carried out in the field of biomarkers for endometriosis, the multicenter prospective ENDOmiRNA Saliva Test study (NCT05244668). Ziwig Endotest® provides information as to the presence or absence of endometriosis to healthcare professionals who prescribe the test.

Ziwig Endotest® has obtained the CE label delivered by a notified body (2797).

# Endometriosis, an all too scarcely recognized disease



Endometriosis affects about 10% of women of childbearing age world wide, i.e. about 190 million women. It potentially affects all women and can appear as early as from adolescence<sup>1</sup>.

Endometriosis is characterized by the presence of fragments of endometrial tissue (a mucous membrane lining the interior of the uterus) outside of the uterine cavity (womb). The presence of these tissue fragments is frequently associated with an inflammatory process<sup>10</sup>. During menstruation, endometriosis lesions can react to variations in hormonal levels, leading to acute inflammation, bleeding, severe pain and other symptoms.

Endometriosis can be completely asymptomatic, and is sometimes diagnosed during a consultation prompted by difficulty in conceiving. In fact, a significant proportion of patients with endometriosis also suffer from infertility<sup>10,11,12</sup>.



# A major impact on personal and professional Wellbeing



Endometriosis can profoundly affect quality of life, interfering with social engagement, physical and sexual health, relationships, work and academic performance, mental health, and overall well-being. Over time, these ongoing challenges may lead to significant barriers in achieving life aspirations, such as completing education, advancing in a career, establishing meaningful relationships, or starting a family. As a result, it can reshape an individual's life path in impactful ways<sup>12</sup>.

Moreover, the persistence of pain can induce a phenomenon of hypersensitization: pain perception thresholds decrease over time leading to the onset of chronic pain<sup>13,14</sup>.

This phenomenon can occur at any stage of endometriosis and chronic pain can persist even after lesions have been treated<sup>15</sup>.

# When to use the salivary diagnostic test?

Ziwig Endotest® paves the way for an early diagnosis of endometriosis and represents a major paradigm shift in clinical practice<sup>3</sup>.

## The key advantages of Ziwig Endotest®:

- Clear Positive/negative result<sup>9</sup>
- Rapid diagnosis<sup>9</sup>
- All forms of endometriosis can be diagnosed<sup>3-8</sup>
- High reliability: Sensitivity > 97%, specificity > 93%<sup>9</sup>
- Simple and non-invasive sample collection<sup>9</sup>







# When to prescribe the saliva test?

The prescription of Ziwig Endotest® is intended for patients between the age of 18 and 43 years with symptoms suggestive of endometriosis<sup>9</sup>.

## Suggestive symptoms of endometriosis<sup>16</sup>:

- Chronic pelvic pain
- +/- dysmenorrhea
- +/- deep dyspareunia
- +/- dysuria/painful urination
- +/- dyschezia/painful defecation
- +/- painful rectal bleeding or haematuria during menstruation
- +/- pain in the shoulder tip
- +/- infertility



Ziwig Endotest® is not intended to replace diagnosis by pelvic medical imaging. Medical imaging is complementary to Ziwig Endotest® and remains essential for lesion description, mapping, and classification.


# Instructions for collecting saliva samples

## 30 minutes before sampling:

- Wash your hands
- Do not use if you are ill
- Do not eat or drink
- Do not smoke
- Do not chew gum
- Do not wear lipstick
- Do not brush your teeth




The quantity of liquid must be sufficient and exceed the fill line (bubbles not included).



1. Spit into the tube until the saliva is above the fill line (excluding bubbles).



2. Close the tube tightly by pressing strongly on the funnel cap until you hear a clear click.



3. Unscrew the funnel from the tube.



4. Close the tube with the small cap



5. Shake the capped tube for 5 seconds.

# Treatment and care

Currently, treatment regimens aim to slow down the progression of the disease, combat symptoms and improve quality of life. Furthermore, an optimized treatment strategy can be implemented in case of associated infertility.

In addition to personalized analgesic and supportive care, hormonal treatments, including estrogen-progestin or progestin-only therapies, are used to regulate or suppress menstruation. These treatments are designed to alleviate symptoms and improve patients' quality of life. As a second-line option, GnRH analogues may be prescribed, to further reduce symptoms. To mitigate side effects like hot flashes, this therapy is typically combined with "add-back therapy." Furthermore, as research continues to advance, new therapies are becoming increasingly available to patients, offering more options for effective management and improved outcomes.

Surgical treatment is considered when medical treatments are insufficient to alleviate symptoms. The objective of surgery is to remove lesions and correct anatomical abnormalities caused by the disease. The procedure is most often performed by laparoscopy.



# Sciences



Ziwig Endotest® uses two cutting-edge technologies: Next Generation Sequencing (NGS) and Artificial Intelligence (AI) to reach a robust diagnosis. In fact, 109 different microRNAs are analyzed in each sample to determine the absence or presence of endometriosis<sup>5,6</sup>.

The development of a non-invasive test for the diagnosis of endometriosis has been the subject of intense medical research for many years. More than 100 biomarkers have been evaluated over the last decades. Among them, a new class of molecules first described in 1993, microRNAs, has emerged as a promising option, supported by a growing number of studies across different medical specialties.

MicroRNAs are small non-coding RNAs. They participate in gene expression: when a microRNA binds to its target, a specific messenger RNA, it primarily blocks protein synthesis, this phenomenon is known as silencing.

MicroRNAs are released into the extracellular environment through various transport structures that shield them from being degraded, thereby enhancing their stability. These circulating microRNAs are present in varying concentrations in different body fluids, including blood, urine, breast milk, tears, and saliva<sup>17</sup>.

In recent years, there has been increasing evidence of the contribution of microRNAs in the pathophysiological mechanisms of endometriosis. A direct link has been made between the dysregulation of certain microRNAs and the development of endometriosis lesions<sup>18</sup>.

# The studies on Ziwig Endotest® have been published in international peer-reviewed scientific journals.

Find all the publications and their abstracts on the website [www.ziwig.com](http://www.ziwig.com)



August 2024

Clues to revising the conventional diagnostic algorithm for endometriosis.

[Read article](#)



June 2023

Validation of the Salivary miRNA Signature of Endometriosis — Interim Data

[Read article](#)



May 2023

Value of non-coding RNAs to assess lymph node status in cervical cancer

[Read article](#)



November 2022

Saliva microRNA signature to diagnose endometriosis: A cost-effectiveness evaluation of the Endotest®

[Read article](#)



October 2022

Endometriosis-associated infertility diagnosis based on saliva microRNA signatures

[Read article](#)



July 2022

A Bioinformatics Approach to MicroRNA-Sequencing Analysis Based on Human Saliva Samples of Patients with Endometriosis

[Read article](#)



May 2022

Endometriosis Associated-miRNome Analysis of Blood Samples: A Prospective Study

[Read article](#)



March 2022

MicroRNome analysis generates a blood-based signature for endometriosis

[Read article](#)



January 2022

Salivary MicroRNA Signature for Diagnosis of Endometriosis

[Read article](#)



January 2022

Clues for Improving the Pathophysiology Knowledge for Endometriosis Using Serum Micro-RNA Expression

[Read article](#)



January 2022

Machine learning algorithms as new screening approach for patients with endometriosis

[Read article](#)



# About Ziwig

Ziwig® is an innovative French biotech specialising in salivary RNA analysis and AI. It operates as an ecosystem at the crossroads of several medical, scientific and digital disciplines. Ziwig® is strongly committed to disruptive innovations that support the transformation of healthcare systems towards more effective, more humane and more accessible precision medicine.

[www.ziwig.com](http://www.ziwig.com)

## Acknowledgements

Ziwig®'s ambition to foster research and treatment of endometriosis and to change the way we understand this debilitating disease that affects so many women around the world, would never have succeeded without all the women and men who believed in our project.

They fought alongside us to make this innovation possible. On behalf of all those who will benefit from our progress, we thank them!



Ziwig Endotest® wins  
the Prix Galien 2024

Ziwig  
19 rue Riboud  
69003 Lyon  
France

ZE\_BR12P\_EN\_Web\_V02\_140125

\*Omnigene Oral is an IVD for the collection and stabilisation of nucleic acids in saliva.  
Manufacturer of the Omnigene Oral : DNA Genotek Inc.

# Bibliographic references

1. Zondervan, Krina T et al. "Endometriosis." The New England journal of medicine vol. 382,13 (2020):1244-1256. doi:10.1056/NEJMra1810764
2. Missmer, Sims OT, Gupta J, Missmer SA, Aninye IO. Stigma and Endometriosis: A Brief Overview and Recommendations to Improve Psychosocial Well-Being and Diagnostic Delay. Int J Environ Res Public Health. 2021 Aug3;18(15):8210. doi: 10.3390/ijerph18158210. PMID: 34360501; PMCID:PMC8346066.
3. Machine learning algorithms as new screening approach for patients with endometriosis. Bendifallah S & al. Sci Rep. 2022 Jan 12;12(1):639. doi: 10.1038/s41598-021-04637-2. [https://www.ncbi.nlm.nih.gov/pmc/articles/PMC8755739/pdf/41598\\_2021\\_Article\\_4637.pdf](https://www.ncbi.nlm.nih.gov/pmc/articles/PMC8755739/pdf/41598_2021_Article_4637.pdf)
4. Clues for Improving the Pathophysiology Knowledge for Endometriosis Using Serum Micro-RNA Expression. Dabi Y & al. Diagnostics (Basel). 2022 Jan 12;12(1):175. doi:10.3390/diagnostics12010175. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC8774370/pdf/diagnostics-12-00175.pdf>
5. Salivary MicroRNA Signature for Diagnosis of Endometriosis. Bendifallah & al. J Clin Med. 2022 Jan 26;11(3):612. doi: 10.3390/jcm11030612. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC8836532/pdf/jcm-11-00612.pdf>
6. Bendifallah, Sofiane et al. "Validation of a Salivary miRNA Signature of Endometriosis — Interim Data" NEJM Evidence. 09 June. 2023. doi: 10.1056/EVIDoa2200282
7. Endometriosis Associated-miRNome Analysis of Blood Samples:A Prospective Study Bendifallah S & al. Diagnostics. 2022; 12(5):1150. <https://doi.org/10.3390/diagnostics12051150>
8. Pascoal E, Wessels JM, Aas-Eng MK, Abrao MS, Condous G, Jurkovic D, Espada M, Exacoustos C, Ferrero S, Guerriero S, Hudelist G, Malzoni M, Reid S, Tang S, Tomassetti C, Singh SS, Van den Bosch T, Leonardi M. Strengths and limitations of diagnostic tools for endometriosis and relevance in diagnostic test accuracy research. Ultrasound Obstet Gynecol.2022 Sep;60(3):309-327. doi: 10.1002/uog.24892. PMID: 35229963.
9. Instructions for use of the IVD device
10. Saunders PTK, Horne AW. Endometriosis: Etiology, pathobiology, and therapeutic prospects. Cell. 2021;184(11):2807-2824. doi:10.1016/j.cell.2021.04.041
11. International working group of AAGL, ESGE, ESHRE and WES et al. "An International Terminology for Endometriosis, 2021."Journal of minimally invasive gynecology vol. 28,11 (2021): 1849-1859. doi:10.1016/j.jmig.2021.08.032
12. Missmer SA, Tu FF, Agarwal SK, et al. Impact of Endometriosis on Life-Course Potential: A Narrative Review. Int J Gen Med. 2021;14:9-25. Published 2021 Jan 7. doi:10.2147/IJGM.S261139
13. As-Sanie, Sawsan et al. "Changes in regional gray matter volume in women with chronic pelvic pain: a voxel-based morphometry study." Pain vol. 153,5 (2012): 1006-1014. doi:10.1016/j.pain.2012.01.032
14. As-Sanie, Sawsan et al. "Functional Connectivity is Associated With Altered Brain Chemistry in Women With Endometriosis-Associated Chronic Pelvic Pain."The journal of pain vol. 17,1 (2016): 1-13. doi:10.1016/j.jpain.2015.09.008
15. Coxon, Lydia et al. "Pathophysiology of endometriosis-associated pain: A review of pelvic and central nervous system mechanisms." Best practice & research. Clinical obstetrics & gynaecology vol. 51 (2018): 53-67. doi:10.1016/j.bpobgyn.2018.01.014
16. Becker, Christian M et al. "ESHRE guideline: endometriosis." Human reproduction open vol. 2022,2 hoac009. 26 Feb. 2022, doi:10.1093/hropen/hoac009
17. Panir K, Schjenken JE, Robertson SA, Hull ML. Non-coding RNAs in endometriosis: a narrative review. Hum Reprod Update. 2018;24(4):497-515. doi:10.1093/humupd/dmy014
18. Dana, Parisa M et al. "Circular RNA as a potential diagnostic and/or therapeutic target for endometriosis." Biomarkers in medicine vol. 14,13 (2020): 1277-1287. doi:10.2217/bmm-2020-0167